



# BRIGHT

BEACON40<sup>®</sup>

Smart Home Light Systems

---

Setup Guide

*Lighting the way to bright memories*

# Setup is simple



Plug in the lamp(s)



Press the power button(s)



If you have more than one lamp  
press and hold the pairing button

## WHAT'S IN THE BOX

### **BEACON40® Personal**

- 1 Base lamp
- 1 Power cord



### **BEACON40® Surround**

- 1 Base lamp
- 1 Satellite lamp
- 2 Power cords





**Use the controls on the Beacon40 Base lamp to set your preferences.**

**1 POWER**

Turn your lights on and off.

**2 TIMER**

Set your lights to turn off automatically after 1, 2 or 3 hours.

**3 TONE**

Set your lights to: warm gold, white or cool blue. You may prefer a cooler tone during the daytime and a warmer tone in the evenings.

**4 BRIGHTNESS**

Turn the knob to your desired luminosity.




## For Beacon40 Surround

### 5 PAIRING

Press and hold the button to pair the Base lamp with the Satellite lamp.

## RECOMMENDED USE

### Choose a Comfortable Location



- Set the light(s) 3-15 feet from your location
- Adjust your preferred brightness and tone
- Ensure the lights are unobscured and within your peripheral vision.

### Create a Simple Routine



- Plan for 1 or more hours of passive light therapy
- Use the lights every day
- OPTIONAL: Set the smart home app to automate your daily routine.

Beacon40 lights are designed to fit comfortably into any home or office. Use the lights while you go about your daily routines like working, watching TV or reading.

# SMART HOME SETUP

---



Download the Tuya Smart app and follow the onscreen instructions for sign up.

---

Devices to be added: 1



BEACON40

Do not add

Go to add

Turn your Beacon40 light on and select Manually Add Device. Select Beacon40 to add it to the application.



Enter the password for the 2.4 GHz Wi-Fi network.



Use voice assistant to set daily routines.

---

You can connect Tuya Smart to Google Home, Amazon Alexa or Apple Home.

For complete smart home setup instructions, please visit [getmebright.com/smart-home](https://getmebright.com/smart-home)

## PRECAUTIONS

Beacon40 is designed for use as a complement to your doctor-prescribed therapy. You should not use this product if you have a history of seizures or epilepsy, or if you have been diagnosed with photosensitive epilepsy. Consult a doctor before use if you believe you may be prone to photosensitive epileptic seizures.

This product is not a medical device. BRIGHT services and products are not intended to be used as a medical device. Furthermore, the services and the products are neither regulated nor approved by the U.S. Food and Drug Administration, and are not designed to detect or prevent causes of Alzheimer's disease or any other medical condition.

BRIGHT services and products are not a substitute for medical care or adult supervision. By using BRIGHT products you acknowledge, understand and agree that your use of the services and the product is entirely at your own risk.

A very small percentage of people may experience a seizure when exposed to flashing lights. Even people who have no history of seizures or epilepsy may have an undiagnosed condition that can cause these "photosensitive epileptic seizures". These seizures may have a variety of symptoms including lightheadedness, altered vision, twitching of the eyes or face, jerking or shaking of arms or legs, disorientation, confusion, or momentary loss of awareness. Seizures may also cause loss of consciousness or convulsions that can lead to injury from falling down or striking nearby objects. If you experience any of these symptoms, turn off the light immediately and consult a doctor.

According to IEC62471-1:206-07 Ed. 1.0 photobiological safety of lamps and lamp systems, sources of optical radiation are classified into risk groups subject to their potential photobiological hazard. Beacon40 is GROUP1 (Low Risk). No photobiological hazard under normal behavioral limitations.

BRIGHT product Beacon40 device risk group is defined in IEC 62471-1:206-07 Ed. 1.0.

Viewing the lights closer than 200mm distance should be avoided. Please follow proper installation and use recommendations to avoid related risks.

For additional information please visit  
**[getmebright.com](http://getmebright.com)** or contact us at  
**[info@getmebright.com](mailto:info@getmebright.com)**